



Fourier transform infrared measurement of the HCl spectrum

Our experiment: rovibrational levels of HCl

The FTIR bench is a versatile platform with several different types of attachments For detection of ATR, IRRAS, transmission and microscopy. We are using a Bruker Tensor 22 FTIR. The second slide shows the sample compartment. The beam passes through the middle of the sample chamber







BRUKER

PIKE
MTEC

TENSOR 27


System Certification
 Validated Non-Validated
Instrument Name: Tensor 27
Program: FTIR
Calibration: None
Date of use: 10/1/2013
PIKE Part No: 1076 420-2839 v01.0712
PIKE System: www.pike-technologies.com

Our experiment: rovibrational levels of HCl

The FTIR bench is a versatile platform with several different types of attachments for detection of ATR, IRRAS, transmission and microscopy.

We use a 10 cm quartz cuvette. Why quartz? What materials are used for FTIR? What do they cost? Do a little web search to answer this question (with any example that you find reasonable).

Transmission sample holder and quartz sample cell



Cuvette

Sample Holder

Our experiment: rovibrational levels of HCl


The FTIR bench is a versatile platform with several different types of attachments for detection of ATR, IRRAS, transmission and microscopy.

We use a 10 cm quartz cuvette. Why quartz? What materials are used for FTIR?

Quartz is much cheaper than other materials such as calcium fluoride or zinc selenide. It absorbs very strongly near 1100 cm^{-1} due to the $-\text{Si-O-Si}-$ modes. Even the windows of the cell absorb 100% of the light below 2200 cm^{-1} . We can use this cuvette only because HCl absorbs at ca. 2900 cm^{-1} .

Add concentrated HCl to the cuvette using a glass pipette. How can we detect HCl?

Prepare HCl sample using concentrated acid



Our experiment: rovibrational levels of HCl

The FTIR bench is a versatile platform with several different types of attachments for detection of ATR, IRRAS, transmission and microscopy.


We use a 10 cm quartz cuvette. Why quartz? What materials are used for FTIR?

Quartz is much cheaper than other materials such as calcium fluoride or zinc selenide. It absorbs very strongly near 1100 cm^{-1} due to the $-\text{Si-O-Si}-$ modes. Even the windows of the cell absorb 100% of the light below 2200 cm^{-1} . We can use this cuvette only because HCl absorbs at ca. 2900 cm^{-1} .


Add concentrated HCl to the cuvette using a glass pipette. How can we detect HCl?

Concentrated HCl has a very high vapor pressure. In a closed cell it vaporizes and fills the cell with enough HCl to easily detect the rovibrational bands.


Background, **silica** and **HCl**




Focus on **silica** and **HCl** to obtain absorption spectrum




Absorbance of HCl sample




Absorbance of HCl sample



Comparison of absorbance of two HCl samples




The raw data obtained from the baseline corrected spectrum shown in the Figure



$$\tilde{\nu}_{obs} = \tilde{\nu}_o + 2(B_e - \alpha)m - \alpha m^2$$

Linear fit of HCl rovibrational spectra



These data were fit using linear least squares. Using the values in wavenumbers we fit the linear function,

$$\tilde{\nu} = a_0 + a_1 m + a_2 m^2$$

Where

$$a_0 = \tilde{\nu}_0$$

And

$$a_1 = 2\tilde{B} - 2\alpha$$

$$a_2 = -\alpha$$

The fit parameters are:

n	0	1
a_n	17.36	2908.00
σ_n	0.3258	1.5515

From these fit parameters we can obtain the Q band wave number $\tilde{\nu}_0$ and the rotational constant \tilde{B} cm⁻¹.

Centrifugal Distortion

Centrifugal distortion arises from the changing bond length of the molecule during rotation. Clearly the rigid rotator approximation is not applicable to a real molecule and the chemical bond acts like a spring that alters the distance between the nuclei. The moment of inertia increases with the bond length squared and therefore the rotational constant decreases as the bond length increases. We can express this phenomenon as a correction to the rotational energy.

$$F(J) = BJ(J + 1) - DJ^2(J + 1)^2$$

where D is the centrifugal distortion constant. D is related to the vibrational wavenumber,

Theoretical basis for quadratic fit

$$D = \frac{4B^3}{\omega^2}$$

When the above factors are accounted for, the actual energy of a rovibrational state is

$$S(v, J) = \nu_0 \left(v + \frac{1}{2} \right) + B_e J(J + 1) - \alpha_e \left(v + \frac{1}{2} \right) J(J + 1) - D_e [J(J + 1)]^2$$

